

CAMMY®

USB INTRAORAL CAMERA



CAMMY Software

CAMMY®

- Directly connects to a computer via the USB port.
- Self-sufficient, CMOS technology, digital image processing.
- Four white LED lamp illumination to provide maximum realistic color reproduction.
- Built-in lens protector designed to prevent dust and scratches to optical components.
- Mini-head design allows for easy access to the third molar.
- Lightweight (only 92 grams), ergonomically designed handpiece.
- User-friendly software. Easy to use and integrate into patient management software.
- Digital images can be easily stored, retrieved and transmitted electronically.
- Disposable sleeves included to prevent cross contamination.

ISO 13485 **GMP**

DENTAMERICA®



USB INTRAORAL CAMERA

1/4" HIGH RESOLUTION COLOR CMOS

LENS PROTECTOR TO PREVENT ACCIDENTAL LENS SCRATCH

4 LED WHITE LIGHT, COLOR TEMPERATURE APPROX. 5000°K

ERGONOMIC COMFORT GRIP OFFERS SMOOTH AND EASY HANDLING

Reorder NO. 336
ZOOM LENS

- Use the zoom lens to magnify and examine caries, root canal access etc.
- Attach the surface viewer ring to examine occlusal surfaces.



MINI-HEAD DESIGN ALLOWS FOR EASY ACCESS TO THE THIRD MOLAR

ERGONOMIC COMFORT GRIP OFFERS SMOOTH AND EASY HANDLING

FREEZE BUTTON

ON/OFF BUTTON

LIGHTWEIGHT HANDPIECE

SPECIALLY DESIGNED HANDPIECE CONTOURED TO INTEGRATE INTO ALL DENTAL UNITS

HIGHLY FLEXIBLE CABLE (NO FRAGILE FIBER OPTIC CABLE) SINGLE PLUG MAKES CAMERA PORTABLE AND MOBILE

SPECIFICATIONS

INTRAORAL USB CAMERA

Interface	USB 2.0
Sensor	1/4" OVT CMOS Sensor
Effective Pixels	0.3Mega Pixels
Max. Resolution	640*480
Resolution Combination	640*480,352*288,320*240, 176*144,160*120
Color Depth	24-bit(16.7 million color)
Max. Frame Rate	30fps
White Balance	Auto
Auto Exposure (AE)	Auto

Motion JPEG Compression	Yes(Smooth Video)
Rating Voltage	5V
Rating Power Consumption	1500mW
Rating Input Current	300mA
USB Cable Length	260 cm
USB Type	USB_A
Weight(g)	92g
Dimensions(L)X(W)X(H)(mm)	182mm*24mm*20mm
Working Temperature	-10°~70°C
Working Humidity	20%~80%

DENTAMERICA®

18688 E. San Jose Ave.
City of Industry, CA 91748 U.S.A.
Phone: (626) 912-1388
Fax: (626) 913-0510
www.dentamerica.com
www.dentamerica.asia

DISTRIBUTOR