

ENDOMAX[®] *BRITE*

Cordless Endodontic Handpiece
with LED Illumination



- Built-in LED illuminates the oral cavity for easy operation.
- Cordless convenience offers complete portability.
- Easy to read LCD display.
- Heavy-duty high torque motor.
- Automatic forward/reverse functions.
- Switch from rotary to reciprocating with a press of a button.
- Five torque settings and nine rotation speeds available for different applications.
- High capacity lithium-ion battery.

ISO 13485




2460

ENDOMAX[®] BRITE


Cordless Endodontic Handpiece with LED Illumination




Auto-Reverse Forward:



Reciprocating Mode:



Auto-Reverse Stop:



SPECIFICATIONS

INPUT	AC 100V~240V, 50/60Hz
FUSE	1A 250V
BATTERY TYPE	3.6V 2040mAh
DIMENSIONS (HANDPIECE)	L201 xW26 xH31 mm
DIMENSIONS (CHARGING STATION)	L105 xW85 xH59 mm
WEIGHT (HANDPIECE)	160g
WEIGHT (CHARGING STATION)	125g

DENTAMERICA[®]

18688 E. San Jose Ave.
 City of Industry, CA 91748 U.S.A.
 Phone: (626) 912-1388
 Fax: (626) 913-0510
 www.dentamerica.com
 www.dentamerica.asia

DISTRIBUTOR