

Steam Sterilizer

Autoclave



Steam Sterilizer

- Easy operation with cycle timer and temperature gauge displaying both chamber temperature and sterilization pressure.
- Stainless steel boiler chamber and sterilization box provide long lasting durability.
- Large volume chamber allows for efficient sterilization.
- Safety thermostat control automatically prevents unit from overheating.
- Double pressure adjustments can auto-vent excessive internal pressure.
- Audible tone followed by auto shutdown when sterilization program is complete.
- Back drain valve design allows for easy cleaning and maintenance.

Steam Sterilizer

Autoclave

C. Door Compel Closing Fixed Axle



A. Exhaust Valve

B. Safety Valve

H. Drain Valve

After replacing the used water, close the drain valve.



D. Inside Pot Door Handle

E. Pressure Indicator



F. Timer



G. Overheating Protection Device
Thermostat control



SPECIFICATIONS

Chamber Capacity	23cm × 41cm (depth)
INSIDE CAPACITY	16.20 (L)
OUTSIDE SIZE	47cm × 33cm × 41cm
STERILIZE TEMPERATURE	Highest 127°C, including under temperature
STERILIZE PRESSURE	Highest 1.5kg/cm ² , including under pressure
STERILIZE TIME	30 minutes
HEATER	800~1,000W
MATERIAL QUALITY	Stainless steel SUS#304
	Boiler chamber surface with burnish
ELECTRIC CURRENT	AC 220V, 7A

DENTAMERICA®

18688 E. San Jose Ave.
City of Industry, CA 91748 U.S.A.
Phone: (626) 912-1388
Fax: (626) 913-0510
www.dentamerica.com
www.dentamerica.asia

DISTRIBUTOR

DENTAMERICA®