

# MICROMAX®

## Cordless Prophylaxis Handpiece



MICROMAX®


- Cordless convenience offers complete portability.
- Heavy-duty high torque brushless motor.
- Forward/reverse functions.
- Reciprocating setting available without use of additional attachment.
- Nine rotation speeds available for different applications.
- Ergonomic design provides comfort and easy operation.
- High capacity lithium-ion battery provides stable and prolonged usage.
- Offers quick and easy coupling.
- Low battery warning light.
- Four different attachments available to choose from.

ISO 13485 QMS CE 2460

DENTAMERICA®

# MICROMAX<sup>®</sup>

## CORDLESS PROPHYLAXIS HANDPIECE



### Optional Attachments



Latch Attachment



Screw Attachment



Snap-on Attachment



Straight Attachment



### SPECIFICATIONS:

INPUT	AC 100V~240V 50/60Hz
FUSE	1A 250V
BATTERY TYPE	3.6V 2040mAh
DIMENSIONS (HANDPIECE)	L192 xW26 xH29 mm
DIMENSIONS (CHARGING STATION )	L105xW85xH59mm
WEIGHT (HANDPIECE)	150g
WEIGHT (CHARGING STATION )	125g

### Order No.

177401 Latch type

177402 Screw type

177403 Snap-on type

177408 Straight type

### DENTAMERICA<sup>®</sup>

18688 E. San Jose Ave.  
City of Industry, CA 91748 U.S.A.  
Phone: (626) 912-1388  
Fax: (626) 913-0510  
www.dentamerica.com  
www.dentamerica.asia

DISTRIBUTOR