

# ROTEX® 780E *PLUS*

Compact Electric E-Type Handpiece Unit



- New control panel with improved ergonomic design.
- Compact size with small foot control.
- Forward/Reverse and Hand/Foot control options.
- E-Type handpiece with straight and contra angle attachments.
- Ball bearing high torque electric motor.
- 0 to 35,000 RPM variable speed.
- Automatic cooling fan for continuous high speed performance.
- Unique compact foot control designed to be activated at any angle and accessible from 360 degrees.
- Application for grinding, polishing, milling of precious metals, ceramics, acrylics, and steel.

**DENTAMERICA®**

ROTEX® 780E *PLUS*

# ROTEX<sup>®</sup> 780E PLUS

## Compact Electric E-Type Handpiece Unit



### 1. Straight & Contra-angle Attachment

Sterilisable up to 135 °C.

Exclusive dental quality with best treatment experience; suit for grinding, polishing, milling of precious metals, ceramics, acrylics and steel.



### 2. Foot Control

Unique foot control designed to be accessible from 360 degrees.



### 3. New Control Panel

Blue LED display makes selecting mode and rotation direction clearly.



Hand Control Mode



Foot Control Mode



Reverse Rotation

Forward Rotation

#### SPECIFICATIONS:

- **INPUT:** AC 110/220V 50/60Hz
- **DIMENSIONS:** L136 X W136X H 75 mm
- **SPEED:** 0-35,000 RPM
- **WEIGHT:** 1.4 kg

#### DENTAMERICA<sup>®</sup>

18688 E. San Jose Ave.  
City of Industry, CA 91748 U.S.A.  
Phone: (626) 912-1388  
Fax: (626) 913-0510  
[www.dentamerica.com](http://www.dentamerica.com)  
[www.dentamerica.asia](http://www.dentamerica.asia)

DISTRIBUTOR

DENTAMERICA<sup>®</sup>